




Rapid Culture Method

Salmonella Precis™

A quick and easy method for the enrichment, detection and confirmation of *Salmonella* species from food, animal feed and environmental samples.

- Validated by AFNOR to ISO 16140 standard
- Simple and easy procedure – no specialised equipment required
- Single 18-hour enrichment
- Single sample transfer
- Single 24-hour plate incubation
- Quick and convenient confirmation: Oxoid Salmonella Latex Test or ISO 6579:2002 standard tests
- Reduced time to result: 2 days compared with up to 5 days for standard culture methods



Introduction

The Oxoid Salmonella Precis™ method combines the benefits of ONE Broth-Salmonella, *Brilliance*™ Salmonella Agar and the Oxoid Salmonella Latex Test to reduce time to result over conventional culture methods.

ONE Broth-Salmonella is a highly nutritious medium for the recovery and growth of salmonellae while inhibiting competing organisms. The growth promoter in the medium allows the recovery of stressed *Salmonella* cells, even when present in very low numbers.

Brilliance Salmonella is the first in a new class of chromogenic media to incorporate novel Inhibigen™ technology. This new technology improves recovery of *Salmonella* by reducing background flora. Chromogens aid easy identification and differentiation by producing brightly coloured colonies.

The Oxoid Salmonella Latex Test provides a quick and easy method for confirmation of *Salmonella* species from culture media.

AFNOR Validation

The Salmonella Precis method has been validated and approved by AFNOR according to ISO 16140 standard against the reference method ISO 6579:2002 standard for the detection of *Salmonella* in food, animal feed and environmental samples.

For flexibility, confirmation was validated using both Oxoid Salmonella Latex Test and the tests outlined in ISO 6579:2002. Alternatively, biochemical panels such as Microbact GNB 24E or RapID ONE Panel, may be used.

AFNOR validation certificate number **UNI 03/06 – 12/07** (available in PDF format from the AFNOR website www.afnor-validation.com).

Reactions on Brilliance Salmonella Agar

| | Colony colour | | |
|--|--|---|--|
| | Purple | Blue | Colourless |
| Enzyme targeted by chromogen | <i>Salmonella</i> (including Lactose positive <i>Salmonella</i>) | <i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i> | <i>Citrobacter</i> , other bacteria and yeasts |
| Esterase | + | -/+ | - |
| β-glucosidase | - | + | - |
| <i>E. coli</i> and other bacteria and yeasts are inhibited by the combination of Inhibigen and other selective agents in the medium. | | | |

| Enrichment | SIZE/FORMAT | ORDER CODE |
|---|-------------|------------|
| ONE Broth-Salmonella in Bottles | 10 x 225ml | BO1096S* |
| ONE Broth-Salmonella in ReadyBags | 3x3 litres | FR60101* |
| ONE Broth-Salmonella Base | 500g | CM1091B |
| ONE Broth-Salmonella Supplement for 225ml | 10 vials | SR0242E |
| ONE Broth-Salmonella Supplement for 2.25 litres | 10 vials | SR0242B |

Plating

| | | |
|---|-----------|----------|
| <i>Brilliance</i> ™ Salmonella (ready to use 90mm plates) | 10 plates | PO5098A* |
| <i>Brilliance</i> ™ Salmonella Agar Base | 500g | CM1092B |
| Salmonella Selective Supplement for 500ml | 10 vials | SR0194E |

Confirmation

| | | |
|--|-------------|-----------|
| Salmonella Latex Test | 100 tests | FT0203A |
| Nutrient Agar | 500g | CM0003B |
| Microbact GNB 24E | 40 tests | MB1131A |
| Microbact GNB 24E | 80 tests | MB1074A |
| RapID ONE Panel | 20 panels | R8311006* |
| Triple Sugar Iron Agar | 500g | CM0277B |
| Urea Agar Base | 500g | CM0053B |
| Urea 40% Solution (for 100ml medium) | 10 x 5ml | SR0020K |
| Lysine Decarboxylation Broth Tablets (each make 5ml) | 100 tablets | CM0308S |
| Spot Indole Reagent (DMACA) | 25ml | R21245* |
| Salmonella O and H Agglutinating Sera various codes - see Product List | | |

Quality Control Organisms

| | | |
|---|----------|--------|
| <i>Salmonella</i> Typhimurium ATCC® 14028™† | 10 loops | C6000L |
| <i>Staphylococcus aureus</i> ATCC® 25923™† | 10 loops | C7010L |
| <i>Klebsiella pneumoniae</i> ATCC® 13883™† | 10 loops | C7037L |
| <i>Enterococcus faecalis</i> ATCC® 29212™† | 10 loops | C7030L |
| <i>Escherichia coli</i> ATCC® 25922™† | 10 loops | C7050L |

* Check code and availability with your local Oxoid Representative

For more information about these products, please visit www.oxoid.com

Protocol for Salmonella Precis Method

Day 0: Enrichment

25g or 25ml of sample
+
225ml ONE Broth-Salmonella



Incubate for
16-24 hours at 42°C

Day 1: Plating


Using a 10µl microbiological
loop inoculate a single *Brilliance*
Salmonella plate




Incubate for
22-26 hours at 37°C

Day 2: Results

If present, select a well isolated purple coloured
colony and test using the Oxoid Salmonella Latex
Test. Alternatively, confirm purple colonies using
standard ISO methods.



Select purple
colonies for confirmation



DEDICATED TO MICROBIOLOGY

Oxoid, Wade Road, Basingstoke,
Hants RG24 8PW

Tel: +44 (0) 1256 841144
Fax: +44 (0) 1256 329728
Email: oxoid.info@thermofisher.com

www.oxoid.com
www.thermofisher.com