



Rapid Culture Method

Salmonella Precis

A quick and easy method for the enrichment, detection and confirmation of *Salmonella* species from food, animal feed and environmental samples.

- Validated by AFNOR to ISO 16140 standard
- Simple and easy procedure – no specialized equipment required
- Single 18-hour enrichment
- Single sample transfer
- Single 24-hour plate incubation
- Quick and convenient confirmation: Oxoid Salmonella Latex Test or ISO 6579:2002 standard tests
- Reduced time to result: 2 days compared with up to 5 days for standard culture methods



Introduction

The Oxoid Salmonella Precis™ method combines the benefits of ONE Broth-Salmonella, *Brilliance*™ Salmonella Agar and the Oxoid Salmonella Latex Test to reduce time to result over conventional culture methods.

ONE Broth-Salmonella is a highly nutritious medium for the recovery and growth of salmonellae while inhibiting competing organisms. The growth promoter in the medium allows the recovery of stressed *Salmonella* cells, even when present in very low numbers.

Brilliance Salmonella Agar is the first in a new class of chromogenic media to incorporate novel Inhibigen™ technology. This new technology improves recovery of *Salmonella* by reducing background flora. Chromogens aid easy identification and differentiation by producing brightly coloured colonies.

The Oxoid Salmonella Latex Test provides a quick and easy method for confirmation of *Salmonella* species from culture media.

Protocol for Salmonella Precis Method


Day 0: Enrichment

25g or 25mL of sample + 225mL ONE Broth-Salmonella
Incubate for 16–20 hours at 42°C




Day 1: Plating

Using a 10µL microbiological loop inoculate a single *Brilliance* Salmonella plate. Incubate for 22–26 hours at 37°C



Day 2: Results

If present, select a well isolated purple coloured colony and test using the Oxoid Salmonella Latex Test. Alternatively, confirm purple colonies using standard ISO methods.



Select purple colonies for confirmation

AFNOR Validation

The Salmonella Precis method has been validated and approved by AFNOR according to ISO 16140 standard against the reference method ISO 6579:2002 standard for the detection of *Salmonella* in food, animal feed and environmental samples.

For flexibility, confirmation was validated using both Oxoid Salmonella Latex Test and the tests outlined in ISO 6579:2002. Alternatively, biochemical panels such as Microbact™ GNB 24E or Remel RapID™ ONE Panel, may be used.

AFNOR validation certificate number **UNI 03/06 – 12/07** (available in PDF format from the AFNOR website www.afnor-validation.com).

Reactions on *Brilliance* Salmonella Agar

	Colony colour		
	Purple	Blue	Colourless
Enzyme targeted by chromogen	<i>Salmonella</i> (including Lactose positive <i>Salmonella</i>)	<i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i>	<i>Citrobacter</i> , other bacteria and yeasts
Esterase	+	-/+	-
β-glucosidase	-	+	-

E. coli and other bacteria and yeasts are inhibited by the combination of Inhibigen and other selective agents in the medium.

Enrichment media	SIZE/FORMAT	ORDER CODE
ONE Broth-Salmonella in Bottles	10x225mL	B01096S*
ONE Broth-Salmonella in ReadyBags	3x3 litres	FR60101*
ONE Broth-Salmonella Base	500g	CM1091B
ONE Broth-Salmonella Supplement for 225mL	10 vials	SR0242E
ONE Broth-Salmonella Supplement for 2.25 litres	10 vials	SR0242B

Plating media		
<i>Brilliance</i> Salmonella (ready to use 90mm plates)	10 plates	P05098A*
<i>Brilliance</i> Salmonella Agar Base	500g	CM1092B
Salmonella Selective Supplement for 500mL	10 vials	SR0194E

Confirmation

Oxoid Salmonella Latex Test	30 tests	FT0203A
Nutrient Agar	500g	CM0003B
Microbact GNB 24E	40 tests	MB1131A
Microbact GNB 24E	80 tests	MB1074A
Remel RapID ONE Panel	20 panels	R8311006*
Triple Sugar Iron Agar	500g	CM0277B
Urea Agar Base	500g	CM0053B
Urea 40% Solution (for 100mL medium)	10x5mL	SR0020K
Lysine Decarboxylation Broth Tablets (each make 5ml)	100 tablets	CM0308S
Spot Indole Reagent (DMACA)	25mL	R21245*
Remel Salmonella O and H Agglutinating Sera various codes - see Product List		

Quality Control Organisms – Culti-Loops™

<i>Salmonella</i> Typhimurium ATCC® 14028™†	5 loops	CL6000*
<i>Staphylococcus aureus</i> ATCC® 25923™†	5 loops	CL7010*
<i>Klebsiella pneumoniae</i> ATCC® 13883™†	5 loops	CL7037*
<i>Enterococcus faecalis</i> ATCC® 29212™†	5 loops	CL7030*
<i>Escherichia coli</i> ATCC® 25922™†	5 loops	CL7050*

† ATCC Licensed Derivative™ The ATCC Licensed Derivative Emblem®, the ATCC Licensed Derivative word mark®, and the ATCC catalog marks are trademarks of ATCC. Oxoid Ltd is licensed to use these trademarks and sell products derived from ATCC® cultures.

* Check code and availability with your local Oxoid Representative

For more information about these products, please visit www.oxoid.com



Oxoid and Remel are specialty microbiology brands of Thermo Fisher Scientific. Our products are available worldwide.

www.oxoid.com
Tel: +44 (0) 1256 841144
oxoid.info@thermofisher.com

© 2010, Thermo Fisher Scientific Inc. All other trademarks are the property of Thermo Fisher Scientific, Inc. and its subsidiaries. Copyright to photos held separately. All rights reserved.

Part of Thermo Fisher Scientific

Folio No 1215/MS/11/10

